

## IMMUNOHISTOCHEMISTRY, ANY 3 MARKERS

| IHC MARKERS | RESULT |
|-------------|--------|
|             |        |
|             |        |
|             |        |

SLIDE NO :

SPECIMEN:

CLINICAL HISTORY:

GROSS:

IMPRESSION:

ADVISED:

INTERPRETATION

| RESULT                          | SCORE |
|---------------------------------|-------|
| Non immunoreactive              | 0     |
| Immunoreactive in 1-25% cells   | 1+    |
| Immunoreactive in 26-50% cells  | 2+    |
| Immunoreactive in 51-75% cells  | 3+    |
| Immunoreactive in 76-100% cells | 4+    |

## Comments

(Link comments with requested markers from George program)

### NOTE:

1. Detection system used is Polymer HRP
2. The impression is based on the material submitted and is not a complete **surgical** pathology report.
3. False negative IHC results due to inadequate fixation of the material sent for evaluation cannot be excluded.

### FIXATION REQUIREMENTS:

- A. The volume of formalin fixative should be atleast 10 times the volume of the specimen.
- B. Decalcification solutions with strong acids should not be used.
- C. Specimens should be immersed in fixative within 1 hour of the biopsy / resection procedure (time of removal & time of immersion to be mentioned)
- D. In all resection (large) specimens, the tumour must be bisected prior to immersion in fixative