

OPIATES SCREEN, URINE
(Lateral Flow Immunoassay)

Positive / Negative

Note:

- A single urine test cannot differentiate casual use from chronic drug abuse.
- Urine drug testing cannot predict extent of impairment, drug dose or exact time of use.
- This is a screening test. All positive results must be confirmed by Gas chromatography / Mass Spectrometry (GC/MS).
- Negative results indicate absence of drug in urine or concentrations below the assay cutoff.

Intended use:

1. Drug abuse treatment programs
2. Pain management clinics.
3. Organ transplantation programs
4. Psychiatric programs

Comments

Opiates (Morphine, Heroin, Codeine) are structurally similar to one another. Morphine is used clinically as a powerful analgesic. Urine specimens test positive for 1 to 3 days following Morphine or Codeine use.