

IMMUNOHISTOCHEMISTRY, TTF-1 (THYROID TRANSCRIPTION FACTOR -1)

IHC MARKER	RESULT
TTF-1 (THYROID TRANSCRIPTION FACTOR-1)	

SLIDE NO:

SPECIMEN:

CLINICAL HISTORY:

GROSS:

IMPRESSION:

ADVISED:

INTERPRETATION

RESULT	SCORE
Non immunoreactive	0
Immunoreactive in 1-25% cells	1+
Immunoreactive in 26-50% cells	2+
Immunoreactive in 51-75% cells	3+
Immunoreactive in 76-100% cells	4+

Comments

TTF-1: This is a nuclear transcription factor necessary for the development of thyroid and pulmonary tissue. It is expressed in all types of Thyroid carcinomas except Anaplastic type. It is also present in most lung carcinomas including Small Cell Carcinoma.

NOTE:

1. Detection system used is Polymer HRP
2. The impression is based on the material submitted and is not a complete surgical pathology report.
3. False negative IHC results due to inadequate fixation of the material sent for evaluation cannot be excluded.

FIXATION REQUIREMENTS

- A. The volume of formalin fixative should be at least 10 times the volume of the specimen.
- B. Decalcification solutions with strong acids should not be used.
- C. Specimens should be immersed in fixative within 1 hour of the biopsy / resection procedure (time of removal & time of immersion to be mentioned)
- D. In all resection (large) specimens, the tumour must be bisected prior to immersion in fixative